

Bioline Salmonella Testkit Selecta


New rapid

AFNOR Approved Method for Detection of Salmonella spp



With the Bioline Salmonella Selecta kit you get an accurate detection of salmonella spp in:

- meat
- milk
- feed
- HACCP
- eggs
- fish
- faeces
- swabs



Rapid: Achieve results after 24 hours of enrichment.

Easy-to-use: Simple procedure. All reagents ready to use.

Reliable: The Bioline Salmonella SELECTA kit detects salmonella at levels as low as 1 cfu/25 g sample.

Approvals: Approved by AFNOR and NordVal.


Attractive price: The Bioline Salmonella Selecta kit is offered at a competitive price..

Safe: No hazardous chemicals.

Getting started: Bioline introduces the method in your laboratory within 2-3 hours.

Method: New rapid enrichment procedure followed by ELISA (Enzyme Linked Immuno-Sorbent Assay).

Time schedule:



bioline

Fredericiavej 414 | Tlf: +45 7662 0709 | E-mail: contact@bioline.dk
 DK-7080 Børkop | Fax: +45 7662 0705 | Website: <http://www.bioline.dk>
 Denmark